Various standard and optional accessories meet more surgery requirements




































OP830 Technical Specifications


- · Table Length: 2020mm
- · Table Width (without side rails): 520mm
- · Table Height (without pads): 680mm ~ 940mm
- · Trendelenburg / Reverse: 30°/30°
- · Lateral Tilt (left/right): 20°/20°
- · Back Place Up/Down: 75°/25°
- · Head Plate Up/Down (detachable): 30°/45°
- · Leg Plate Up/Down (detachable): 5°/90°
- · Leg Plate Split (detachable): 180°
- · Kidney-bridge Elevator Up: 136mm
- · Longitudinal Shift (optional): 300mm
- · Weight Capacity: 297kg

Exclusively Imported and Distributed by:




105 Maryland St., Cubao, Quezon City Tel. No.: (02) 726-8715 / (02) 727-7822

E-mail: info@aaapharma.net Website: www.aaapharma.net



Aegistab OP830 **Electric Operating Table** Aegistab OP830 **Electric Operating Table**

Aegistab OP830 offers high load capacity, flexibility & comfort to patients in thoracic surgery, abdominal surgery, neurosurgery, ENT surgery, gynecological & obstetric surgery, and urological surgery.





- · Integrated rechargeable battery supports OP830 operation for 15 days (typical) after fully charged
- Use of AC power doable in case of low battery power, improved emergency response capacity

Built-in electro-hydraulic system, load capacity up to 297 kg, high safety and stability



- · Durable mattress filled with high density foam is easy to clean the cover, simple to detach, comfortable to use, and supports prevention of decubitus
- · Built-in mesa mattress MDF with free X-ray access, conducive to the C-arm film



- · Pedestal lock/unlock control performed by a hydraulic device which makes OP830 easy to move and stable when locked
- Base locked automatically if there is any repositioning, avoiding unintended operation



Handset controller brings convenience when repositioning in various surgeries



Built-in kidney-bridge can be elevated to up to 136 mm and provides displays requested by gall bladder surgery and kidney surgery



Folding upward & downward of head boards & leg boards by cylinder control for accurate positioning



- \cdot T-shape base design ensure higher stability
- · Space pre-reserved for surgeons to work closer to the operating table

